

# **Simultaneous Measurements of Carbon, Hydrogen, Nitrogen, Sulfur, and Oxygen with Thermal/Optical Analysis**

John G. Watson<sup>1</sup>, Gustavo M. Riggio<sup>1</sup>, Xiaoliang Wang<sup>1</sup>, L.-W. Antony Chen<sup>1,2</sup>, Xufei Yang<sup>3</sup>, Jana Diab<sup>4</sup>, Ralf Zimmermann<sup>4</sup>, and Judith C. Chow<sup>1</sup>

<sup>1</sup>Desert Research Institute, Reno, NV

<sup>2</sup>University of Nevada, Las Vegas, NV

<sup>3</sup>Montana Tech, Butte, MT

<sup>4</sup>University of Rostock, Rostock, German

Presented at:

11<sup>th</sup> International Conference on Carbonaceous Particles in the Atmosphere  
Berkeley, CA  
August 11, 2015


# Motivation

- More than 100,000 thermal/optical analyses (TOA) are performed worldwide on quartz-fiber filters each year, including long-term trends networks in the U.S., Canada, and China
- It is desirable to obtain more information from these analyses beyond simple carbon fractions (e.g., organic and elemental carbon [OC and EC]) at no added cost
- New developments in detector technology show potential for expanding the components quantified by thermal methods

# Objectives

- Identify approaches for expanding thermal analysis from carbon (C) to hydrogen (H), nitrogen (N), sulfur (S), and oxygen (O) and their associated compounds
- Demonstrate that the long-term OC/EC trends record can be maintained by detector modifications

# Approach 1: Emulate the AMS for filters (TOA-QMS\*)




\*TOA-QMS: Thermal/Optical Analyzer with Quadrupole Mass Spectrometry

Riggio, G.M. (2015). Development and application of thermal/optical- quadrupole TOA-QMS mass spectrometry for quantitative analysis of major particulate matter constituents., M.S. Thesis, University of Nevada Reno, Reno, NV.

# The IMPROVE temperature program is simplified for testing


## Temperature Steps (in helium atmosphere):

- 80 °C – Desorption of H<sub>2</sub>O
- 580 °C – Combustion of most species (OC4 of IMPROVE\_A)
- 840 °C – Combustion of remaining species (EC3 of IMPROVE\_A)




# Air infiltration during sample insertion is reduced with higher pressure and an argon sheath


Increase sample oven pressure from 5 to 10 psi



Reduce O<sub>2</sub> from 15.5 to 8.5 ppm



Argon sheath chamber





# TOA-QMS spectra are similar , but not identical, to AMS spectra



- Potential causes of differences

- Heating rate

- Particle collection medium

- Thermal desorption

- Ionization








## AMS Spectra from:

Allan et al. (2004). A generalised method for the extraction of chemically resolved mass spectra from aerodyne aerosol mass spectrometer data. *J. Aerosol Sci.*, **35**(7):909-922.

Jimenez et al. (2003). Ambient aerosol sampling using the Aerodyne aerosol mass spectrometer. *J. Geophys. Res.*, **108**(D7):doi:10.1029/2001JD001213.





Signal/response is determined by analysis of quartz filter samples of nebulized  $\text{NH}_4\text{NO}_3$ ,  $(\text{NH}_4)_2\text{SO}_4$ , and oxalic acid solutions.






# Application to ambient samples shows good correlation, but systematic biases

(58 Fresno, CA, samples; Dec 2000 – Feb 2001)




# Concentration variations are similar to those obtained from speciation analyses

(Fresno supersite, 58 samples; Dec 2000 – Feb 2001)







Largest concentration of species between 16:00 and 24:00 Local Standard Time

# Approach 2: Oxidize thermally-evolved products to simpler compounds (TOA-O-QMS/NDIR)








# Existing thermal/optical protocols can be adapted to quantify C, H, N, S, and O



Fresno and Baltimore ambient samples (N=87)


# Instrument signals are linear with C, H, N, S, and O quantities for model compounds



## Calibration compounds:

- Ammonium nitrate:  $\text{NH}_4\text{NO}_3$ ;
- Ammonium sulfate:  $(\text{NH}_4)_2\text{SO}_4$ ;
- Benzoic acid:  $\text{C}_7\text{H}_6\text{O}_2$ ;
- Carbon dioxide:  $\text{CO}_2$ ;
- L-Cystine:  $\text{C}_6\text{H}_{12}\text{N}_2\text{O}_4\text{S}_2$ ;
- Methane:  $\text{CH}_4$ ;
- Nonadecanol:  $\text{C}_{19}\text{H}_{40}\text{O}$ ;
- Pentadecanoic acid:  $\text{C}_{15}\text{H}_{30}\text{O}_2$ ;
- Sulfanilamide:  $\text{C}_6\text{H}_8\text{N}_2\text{O}_2\text{S}$

# NDIR signal is linear with O quantities in calibrated chemicals



## Challenges:

- $H_2O$  bound to filters and particles
- Intrusion of ambient  $O_2$  into the analyzer

Benzoic acid:  $C_7H_6O_2$ ;


Nonadecanol:  $C_{19}H_{40}O$

Pentadecanoic acid:  $C_{15}H_{30}O_2$




# CHNSO concentrations are comparable with other methods


**N (DRI vs. Thermo CHNSO Analyzers)**




**N (DRI CHNSO Analyzer vs. IC)**




**TC (DRI vs. Thermo CHNSO Analyzer)**




**S (DRI CHNSO Analyzer vs. IC)**



**H (DRI vs. Thermo CHNSO Analyzers)**






**O (DRI vs. Thermo CHNSO Analyzer)**



- **Thermo Flash EA1112 CHNS/O Analyzer**
- **Dionex Model ICS-3000 Ion Chromatographs (IC)**

**Fresno and Baltimore ambient samples (N=87)**

# Composition varies between summer and winter (Fresno, California)




Seasonal variability in the CHNS-O composition of Fresno ambient samples (N = 67)

Abundant  $(\text{NH}_4)_2\text{SO}_4$  in summer (Decompose at 200–400 °C; OC2 at 280 °C in 100% Helium) (EC/TC=0.23)

Abundant  $\text{NH}_4\text{NO}_3$  in winter (Dissociation starts at room temperature; OC1 at 140 °C in 100% Helium) (EC/TC=0.28)

# Source profiles vary between smoldering and flaming wood smoke for thermal carbon fractions

Smoldering Wood Smoke (BIOQCF101)




Smoldering wood smoke shows lower EC:TC and higher O:C ratios than flaming smoke.

$$\text{EC}/\text{TC}=0.02$$

Fraction	Molar Ratios	
	H:C	O:C
OC	1.55	0.41
EC	0.93	NA
TC	1.54	0.40


Flaming Wood Smoke (SDKQ073)





$$\text{EC}/\text{TC}=0.25$$

Fraction	Molar Ratios	
	H:C	O:C
OC	1.33	0.27
EC	0.14	NA
TC	1.04	0.20

# Elemental analysis of organic standards shows consistent O/C and H/C relationships




## H/C relationships






Lower N/C ratio may be due to inadequate number of samples and unaccounted N species

# Approach 3: Detect thermal output with photon ionization time-of-flight mass spectrometry (TOA-PI-TOFMS)



Diab et al. (2015). Hyphenation of a EC/OC thermal-optical carbon analyzer to photo ionization time-of-flight mass spectrometry: A new off-line aerosol mass spectrometric approach for characterization of primary and secondary particulate matter. *Atmos. Meas. Tech. Discuss.*, (8):269-308.

Grabowsky et al. (2011). Hyphenation of a carbon analyzer to photo-ionization mass spectrometry to unravel the organic composition of particulate matter on a molecular level. *Anal. Bioanal. Chem.*, **401**(10):3153-3164.




**Soft ionization  
doesn't fragment  
components, mass  
spectra are more  
complex, but  
individual  
compounds are  
quantified**

\*OC1-OC3 are OC fractions evolved at 140, 280, and 480°C in helium atmosphere following IMPROVE\_A protocol



# Distinct temperature/ion profiles are discernable, even without identifying individual compounds




Gasoline Exhaust

Diesel Exhaust

Grabowsky et al. (2011). Hyphenation of a carbon analyzer to photo-ionization mass spectrometry to unravel the organic composition of particulate matter on a molecular level. *Anal. Bioanal. Chem.*, **401**(10):3153-3164.

**These approaches seem to be feasible, but not yet practical**



## **Approaches 1 and 2**

## **Approach 3**



# Mini mass spectrometers are demonstrating sufficient sensitivity for ambient concentrations



Torion Technologies Inc.  
American Fork, UT,  
<http://torion.com/home.html>.



Microsaic Systems. Abingdon, UK,  
<http://www.microsaic.com/home/>



Aston Labs, Purdue University, Lafayette, IN,  
<http://aston.chem.purdue.edu/research/instrumentation/miniature-mass-spectrometers>.

# Challenges for Enhanced Chemical Characterization of Filter Samples

- Perfecting, evaluating, and making more efficient procedures for additional characterization
- Modifying instrumentation and procedures to incorporate more specific analyses methods into long-term chemical speciation networks to obtain more information from existing samples
- Maintaining continuity and consistency with the long-term trends data sets
- Developing more detailed source profiles with these methods for speciated inventories and source apportionment

# Acknowledgements

- U.S. National Science Foundation  
(CHE 1214163)
- National Park Service IMPROVE Carbon  
Analysis Project (C2350000894)